


trident



SOLE

Class B Autoclave 17 and 22 litres



trident



CLASSIFICATION
Medical Device Class IIb
Class B Autoclave fully
conforms to EN13060
standard requirements



DIMENSIONS
46 X 59 X h45 cm



WEIGHT
60 Kg



CYCLES RECORDING
USB connection to extract and
store the sterilization cycles
information in a .txt file



CHAMBER CAPACITY
17 & 22 liters



WATER CONSUMPTION
0.5 liter/cycle



CYCLES
12 cycles including two
customizable cycles




Automatic opening door system.
Front water loading and unloading system.
"On" switch on the front.

The stainless steel chamber is a single block,
without welded seams liable to cause rusting,
leaks or cracks where microorganisms can collect.

Class B autoclave

SOLE

17 & 22 liter



With a 98% air extraction Sole creates the perfect atmosphere to generate entirely water-free dry saturated steam

Because sterilization is not an option

Behind every
great product,
there is a great
brand

trident

Trident S.r.l.
Via Artigiani 4
Castenedolo, 25014 BS
Italy
Tel 39 030 2732485
trident-dental.com