


TECHNICAL CHARACTERISTICS SCANNER	
Supply Voltage	100-240V AC, 50-60 Hz
Input Current	0.1 A (240V)/0.24 A (100 V)
Absorbed Power	24 VA
Device voltage	24 V DC
Device current max	1 A
Dynamic Range	16 bit
Resolution	14 lp/mm
Supply voltage	5 V DC from USB
Connection to PC	USB 2 (A type) cable 2 mt
Maximum absorbed current	275 mA
Frame rate	0.7 fps

TECHNICAL CHARACTERISTICS PHOSPHOR PLATES			
Phosphor Composition	BaSrFBr:Eu		
Luminescence	400 nm		
Phosphor Plate dimensions	Size 0		22x31
	Size 1		24x40
	Size 2		31x41
	Size 3		27x54
Pixel shape	squared		
Pixel size	35 µm		
Image dimensions	Size 0		726x1155
	Size 1		792x1321
	Size 2		1024x1352
	Size 3		891x1783



Trident S.r.l.

Via Artigiani, 4 25014 Castenedolo

Phone +39 030 2732485

info@trident-dental.com

www.trident-dental.com