


## APPLICATIONS

FIELD	INDICATIONS
Prosthodontics	Resin inlays/onlays
	Zirconia copings
	Single-tooth restorations in lithium disilicate, zirconia, all ceramic
	Frameworks and fixed partial dentures in zirconia
	Single implant crowns
	Implant bridges
	Implant supported bars
	Posts and cores
	Partial removable dentures
	Digital smile design
Implantology	Obturator
	Guided implant surgery
Orthodontics	Diagnosis and treatment planning
	Aligners
	Custom made devices

PARAMETER	UNIT	VALUE
Scanning technology	-	Staring scanner
Full jaw scanning time	min	<2
Accuracy	µm	10
Field of View	mm	18x16
Depth of field	mm	0 to 15
Data format	-	STL, PLY, PTY
Dimensions (LxWxH)	mm	216x40x36
Weight	gr	228
Tip size (LxWxH)	mm	83x19x15
Rated voltage	V	AC 100-240
Rated current	A	1.5
Rated Frequency	Hz	50-60

COMPONENT	REQUIREMENTS
Operating System	Windows 7* 64 bit, Windows10* 64 bit
CPU	i7-7700k (3.6GHz; quad-core; 8-thread) or above
RAM Memory	DDR4 2400 16G or above
Hard disk	SSD 240G or above
Graphic card	GTX 1060 , GTX 1660 Ti, RTX 2060 or RTX 3060 (Memory 6 Gb or above)
Graphic drivers version	NVIDIA* Driver 512 - WHQL Version or above
USB Port	3.0 type A

**Trident S.r.l.**

Via Artigiani, 4 25014 Castenedolo

Phone +39 030 2732485

info@trident-dental.com

www.trident-dental.com